



**SURGI-TEC**   
THE ART OF MOVING BONE

**Surgi-Tec offers you a complete range for Trans Palatal Distraction**



### **Trans Palatal Distractor TPD Neo**

Our newest distractor in the Surgi-Tec TPD family

- TPD Neo : the 'slim line' edition
- "Ultra thin" model with all the benefits of our TPD All - in - one
- Unique blocking screw concept

### **Trans Palatal Distractor TPD All - in - one**

- Placement fast and easy
- Mobile joint of abutment plates to distraction unit
- Fixation with "self-drilling" screws, eliminates pre-drilling
- Blocking screw prevents to unscrew the distraction module



### **Trans Palatal Distractor TPD Classic**

- Interchangeable modules, enable large expansion
- Fixation with "self-drilling" screws, eliminates pre-drilling
- Blocking screw prevents to unscrew the distraction module

